© CTS | Ferroperm This document and the information in it is proprietary and is the property of CTS | Ferroperm. It may not be copied or disclosed to a third party or used for any purpose other than that for which it is supplied without the express written consent of CTS | Ferroperm.

Information contained in this document is subject to certain Export Control regulations, specifically of the (choose as appropriate) United Kingdom, European Union, other national jurisdiction and / or the United States International Traffic in Arms Regulations and / or Export Administration Regulations. Each recipient of this document is responsible for ensuring that transfer or use of any information contained in this document complies with all relevant Export Control Regulations.

Piezoelectric single- and multi-element ultrasonic transducers for medical applications

T. Zawada, L.M. Bierregaard, E. Ringgaard, R. Xu, M. Guizzetti





39th International Microelectronics and Packaging Conference September 20-23, 2015 Gdansk, Poland

Outline

1 **Company introduction** 2 **Ultrasonic imaging** 2 **Single element transducers Transducer arrays 3D imaging Conclusions**

Company introduction







Overview

» CTS | Ferroperm is an international global engineering group

 Extreme environment components and subsystems for aerospace, defence and energy markets

» Strong expertise in extreme environment engineering

- High technology products and systems
- Content on virtually every aircraft 63,000 installed base
- » Well-balanced portfolio
 - Commercial/military/other (48%/34%/18%)
 - Original equipment/aftermarket (41%/59%)
- » Global footprint
- » Annual sales of £1.554 billion in 2014
- » FTSE100 company



Original equipment 56% Aftermarket 44%

Global presence

~10 800 employees worldwide



CTS | Ferroperm Denmark

- » CTS | Ferroperm is a manufacturer of piezoelectric materials, components, devices of Ferroperm[™] Piezoceramics and InSensor[™] brands
- » 3-4 million units produced annually
- » Major markets
 - Medical ultrasound
 - Underwater acoustics
 - Acceleration sensors
 - Flow meters
 - NDT
 - Energy Harvesting



2 Medical ultrasonic imaging



Global imaging modalities 2011 vs. 2000



Piezoelectric single- and multi-element ultrasonic transducers for medical applications

Medical imaging and enabling technologies



Ultrasonic imaging – basic idea

- > Ultrasonic imaging is highly detailed and geometrically correct (average speed of sound in a soft tissue is 1540 m/s and varies only 3%)
- The imaging is enabled by differences of acoustic impedance defined as a product of density and speed of sound (Z=pc)
- » The amplitude reflection factor acoustic plane wave incident at the interface of two tissues (Z_1 and Z_2) is then given by $RF=(Z_1 - Z_2)/(Z_1 + Z_2)$



Physical processes in an ultrasonic system



Ultrasonic imaging – basic properties

- » Resolution of ultrasound imaging depends on:
 - Aperture size
 - Center frequency
 - Bandwidth
 - Focal depth
- » Typical frequencies range from 1-15 MHz, corresponding to lateral resolution of 3 mm to 0.2 mm
- » Typically, the best resolution is obtained at the focal depth (due to diffraction)
- » Resolution is also affected by attenuation in media, which also depends on frequency (the higher the frequency, the higher the attenuation)

3A Single element transducers







Typical structure



Modelling of a transducer – KLM model



Matching – optimizing energy transfer

- » The transducer has maximum broadband sensitivity only when properly matched to the driving circuit (electrical matching) and the patient (acoustical matching)
- Typically a passive circuit network is used for electrical matching (targeting 50 Ohm matching)
- » Acoustic matching is typically obtained by using additional matching layers on the front of the transducer
- » Typically $\lambda/4$ layers are used with acoustic impedance between the Z_t = 34 MRayls (typical PZT) and 1.5 MRayls (human body)





Electrical matching

Focusing

- » Focusing is used to increase the spatial resolution in a selected region
- » There are similarities between optics and acoustics, e.g. focal region, focal length
- » However, apertures do not have to be circularly symmetric, lenses can be made out of materials having equivalent refraction index less than one
- » Acoustic focusing can also be achieved by electronic means

After: T. L. Szabo (2014): Diagnostic ultrasound imaging: Inside out *Amsterdam, The Netherlands, Elsevier.*



Piezoelectric single- and multi-element ultrasonic transducers for medical applications

Key

R

n Cı

Cw

Pulse-echo test set-up (round-trip testing)

- » A pulse-echo transducer uses the same element to send out the signal and monitor the returning sound echoes
- The transducer rests on a platform with adjustable angle with respect to the target surface (in order to make the face of the transducer parallel to the face of the steel target)
- » Typically a signal amplitude of 100 V was used, and the period of the pulse was varied to maximise the return signal (~ pulse period corresponding to one half the period of the transducer's resonance frequency)





Drive voltage vs. time

Characteristics of a pulse-echo signal (example)

- The pulse is a high-voltage drive; amplification of the actual signal is turned on after the pulse event is completed
- The main-bang ringdown (MBRD) is a period where high acoustic energy reverberates in the piezoelement
- The signal has a backing echo starting at about 20 µs and a secondary one at about 85 µs
- » The actual target echo arrives at about 60 µs



Frequency domain pulse echo properties



3B PZT Thick film imaging transducers







PZT thick films – InSensor™

» Technology of piezoelectric thick films – enabling deposition and integration of piezoelectric layers (10 μm to 100 μm in thickness) with high lateral resolution (100 μm x 100 μm)



- » Key features of InSensor™ technology
 - Capable of manufacturing miniaturized devices
 - Low prototyping costs
 - High volume production
 - High lateral resolution
 - High frequency
 - High response
 - Piezoelectric material can be deposited on a number of different substrates (compatible with MEMS)

PZT thick film compatibility



Thick-film ultrasonic transducers - concept

- » The characteristics of the printed thick film make it a perfect candidate for medical imaging due to the following:
 - More than 20 % porosity
 - Low acoustic impedance
 - Low permittivity
 - High frequency (more than 20 MHz) easily obtained without machining
- When the thick film is combined with the optimised substrate, a fractional bandwidth above 100 % is measured



Typical structure of a high-frequency acoustic transducer based on thick film

Fabrication method – pad printing



High frequency imaging transducer



P. Maréchal, F. Levassort, J. Holc, L.-P. Tran-Huu-Hue, M. Kosec & M. Lethiecq (2006): High-frequency transducers based on integrated piezoelectric thick films for medical imaging. *IEEE Trans. UFFC*, **53** [8]. cf. F. Levassort, L. P. Tran-Huu-Hue, J. Holc, T. Bove, M. Kosec & M. Lethiecq (2001): High-performance piezoceramic films on substrates for high-frequency imaging. *In* Proc. IEEE Ultrason. Symp. 2011, 1035-1038.

High frequency transducer manufacturing



¹*F*: measured focal distance; DOF_6 : depth of field at -6 dB; Δz_6 : axial resolution, BW_6 : relative bandwidth at -6dB, BW_{20} : relative bandwidth at -20 dB.

F. Levassort, E. Filoux, M. Lethiecq, R. Lou-Møller, E. Ringgaard & A. Nowicki (2006): Curved piezoelectric thick films for high-resolution medical imaging. *In* Proc. IEEE Ultrason. Symp. 2006, 2397- 2400.



High-frequency ultrasonic transducers – acoustic field



Acoustic field of thick film transducers measured by means of Schlieren setup, a) flat transducer, b) focused transducer (in collab. with A. Nowicki, IPPT).

High-frequency ultrasonic transducers – medical imaging



Source: Ulthera video http://www.youtube.com/watch?v=XhafnO0uB_k

4A Multi-element transducers







Why arrays?

- » No moving parts, no mechanical scanning
- » Arrays:
 - Enable electronically controlled beam steering
 - Enable dynamic focusing on receive
 - Improve lateral resolution by adjustment of the length as well as apodization (compensation for finite size of an element)
- » Enable flow imaging



1D array orientation



30 MHz 1D array - an example



Beam steering



Beam focusing



Piezoelectric single- and multi-element ultrasonic transducers for medical applications

4B Thick film arrays for medical imaging







Deposition – screen printing



Concept of thick film based array

- » Fabrication be means of screen printing
- » Special characteristics of PZT thick films open for
 - High frequency > 10 MHz
 - An integrated backing layer
 - Large bandwidth
- » Fully printed design
 - » Huge implication for cost reduction
 - » Arbitrary patterning possible



32 element kerfless thick film ultrasonic array – fully printed



Very good distribution of capacitane over the elements



Piezoelectric single- and multi-element ultrasonic transducers for medical applications

Pulse-echo measurements continued



Bandwidth of each element for two 32-element thick film ultrasonic transducers



Center frequency of each element for two 32-element thick film ultrasonic transducers

- Thick film ultrasonic arrays offer a unique combination of high frequency, good sensitivity and a bandwidth comparable to conventional bulk transducers
- » 32-element transducers generally show acceptable uniformity

Pulse-echo measurements - Crosstalk



- » Crosstalk is measured at time zero by comparing the voltage amplitude for the four elements closest to the exited element
- » The average crosstalk is found to be approx. -38 dB
 - An acceptable level for an ultrasonic transducer

5 3D ultrasonic imaging







3D ultrasonic imaging



Fetus 34 weeks old with bilateral cleft lip visualized using (a) conventional 2-D US imaging and (b) 3-D US volume rendering

Karadayi, K., R. Managuli, and Y. Kim (2009) "Three-Dimensional Ultrasound: From Acquisition to Visualization and From Algorithms to Systems". *IEEE Trans. UFFC, vol. 56, no. 4, April 2009*

3D ultrasonic imaging – scanning strategies



Row-column concept



Acoustic field at transmission and reception



Conclusions







Conclusions

- » Ultrasonic imaging is a versatile, non-invasive and low cost alternative to other diagnostic methods
- » Ultrasound is used mostly for imaging of soft tissue
- Image quality (resolution, contrast, etc.) strongly depends on properties of the ultrasonic transducers
- » Historically, mechanically scanned single element transducers have been used, but these are now being replaced by multi-element array systems
- >> The current development of ultrasonic transducers is aiming at:
 - Higher quality imaging at lower cost
 - High frequency imaging with deeper penetration
 - Integrated probes reducing the impact of expensive and heavy cabling
 - New modalities, e.g. elastography (measurement of elastic properties of tissue)
 - New technologies, e.g. cMUT (capacitive micromachined ultrasonic transducer)
- The thick film process allows for cost-effective manufacturing of linear arrays, and the performance of each element is comparable to the single-element thick films

Acknowledgements

» Danish National Advanced Technology Foundation, contract no. 82-2012-4, FutureSonic platform

- » Dr. Frans Lautzenhiser, Ferroperm, Indiana, USA
- » Rasmus Lou-Møller, Ferroperm, Denmark
- » Dr. Franck Levassort and Prof. Marc Lethiecq, GREMAN, Univ. François-Rabelais de Tours, France
- » Prof. Andrzej Nowicki, IPPT, Warsaw, Poland



The information contained in this document is the property of CTS | Ferroperm and is proprietary and/or copyright material. This information and this document may not be used or disclosed without the express authorization of CTS | Ferroperm. Any unauthorized use or disclosure may be unlawful.

The information contained in this document may be subject to the provisions of the trade compliance regulations (including those regulations governing transfer to a dual national or third country national, export and re-export) of various countries; see the first page for specific requirements. The recipient acknowledges that licences from the applicable regulatory agency or agencies may be required before the recipient may further disclose such information to others, and that such licences may impose restrictions on further disclosure of such information. The recipient agrees to comply with all applicable governmental regulations as they relate to the transfer, export and re-export of information disclosed in this document.